


Specification

DTE BioChemica

A1102

Physical Description:	Solid
Product Code:	A1102
Product Name:	DTE BioChemica
Specifications:	Assay (iodometr.): min. 99 % Solubility (5 %; H ₂ O): clear, colorless pH (0.1 M; H ₂ O; 20°C): 4.0 - 6.0
Hazard pictograms	
WGK:	1
Storage:	2-8°C
Signal Word:	Attention
GHS Symbols:	GHS07
H Phrases:	H315 H319 H335
P Phrases:	P261 P280 P304+P340 P305+P351+P338 P405 P501
Molecular Formula:	C ₄ H ₁₀ O ₂ S ₂
M:	154.25 g/mol

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Specification

DTE BioChemica

A1102

CAS:	6892-68-8
EINECS:	229-998-8
CS:	29309098
<p>Comment</p> <p>Dithioerythritol (DTE) is an isomere of dithiothreitol (DTT). In general, DTE and DTT are exchangeable. DTE has pK values of 9.0 and 9.9 and DTT of 8.3 and 9.5, respectively. Both reagents function in the ionized state only. Therefore the reducing activity is highest in the pH range of 7.0 - 9.5. At lower pH values DTT is more active than DTE. The reducing activity is stopped by acidifying the pH (e. g. acetic acid at pH 3; ref. 3). Proteins may oxidize during gel electrophoresis. This results in diffuse bands (decreased resolution) and a non-reproducible mobility (wrong estimation of molecular size). Reduction of the protein by DTE (0.25 M stock solution in water; 5 mM final concentration) and alkylation by iodoacetamide (0.25 M stock solution; 1 \: 20 diluted in the final concentration) solves the problem (2).</p>	
<p>Bibliography</p> <p>(1)Cleland, W.W. (1964) <i>Biochemistry</i> 3, 480-482Dithiothreitol, a new protective reagent for SH groups. (2)Lane, L.C. (1978) <i>Anal. Biochem.</i> 86, 655-664A simple method for stabilizing protein-sulfhydryl groups during SDS-gel electrophoresis. (3)Jocelyn, P.C. (1987) <i>Methods Enzymol.</i> 143, 246-256Chemical reduction of disulfides.</p>	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340