

Specification

Coomassie® Brilliant Blue R-250 (C.I. 42660)

A1092

Physical Description:	Solid
Product Code:	A1092
Product Name:	Coomassie® Brilliant Blue R-250 (C.I. 42660)
Headline Comment:	® registered trademark of Imperial Industries PLC
Specifications:	E 1 %, 1 cm, λ_{\max} : >300 (pH 7.0) λ_{\max} . (buffer pH 7.0): 554 - 563 nm
WGK:	1
Storage:	RT
Molecular Formula:	$C_{45}H_{44}N_3NaO_7S_2$
M:	825.98 g/mol
CAS:	6104-59-2
EINECS:	228-060-5
CS:	32041200
Comment	<p>Coomassie® Brilliant Blue R-250 is one of the most commonly used stains for proteins, after their separation by polyacrylamide gel electrophoresis. The protein-dye complex has an absorption maximum at 549 nm, the dye without protein at 555 nm (in 0.01 M citrate buffer, pH 3). The intensity in staining of proteins probably depends on the basicity of a protein (5). Per positively charged amino acid approximately 1.5 - 3 molecules of Coomassie® will be bound. This variation complicate the exact protein determination with albumin as a standard, since this protein contains more basic amino acids than many other proteins (5). There do exist many protocols for sensitive staining procedures with Coomassie® (e. g. ref. 3, 4). The sensitivity reaches a limit at 25 ng protein (4). We recommend the following protocol:</p> <p>I. Staining solution: 0.1 % Coomassie® Brilliant Blue R-250 (Prod.-No. A1092) 20 % methanol (or ethanol) 10 % acetic acid</p> <p>The SDS gel (without 'stacking gel') is stained for 1 hour at 60°C or for 2 hours at 50°C or over night at RT.</p> <p>II. Destaining solution: 20 % methanol (or ethanol) 10 % acetic acid</p> <p>Destain the gel for 3 - 4 hours at 50 - 60°C. Add some sponges. Subsequently wash the gel for 15 minuts in water and dry under vacuum at 60°C for 2 - 3 hours.</p>

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Bibliography

(1)Fazekas De St. Groth, S. *et al.* (1963) *Biochim. Biophys. Acta* **71**, 377-391Two new staining procedures for quantitative estimation of proteins on electrophoresis strips. (2)Chrambach, A. *et al.* (1967) *Anal. Biochem.* **20**, 150-154A procedure for rapid and sensitive staining of protein fractionated by polyacrylamide gel electrophoresis. (3)Neuhoff, V. *et al.* (1988) *Electrophoresis* **9**, 255-262Improved staining of proteins in polyacrylamide gels including isoelectric focusing gels with clear background at nanogram sensitivity using Coomassie® Brilliant Blue G-250 and R-250. (4)Choi, J.-K. *et al.* (1996) *Anal. Biochem.* **236**, 82-84Modified Coomassie® Blue staining of proteins in polyacrylamide gels with Bismark brown. (5)Tal, M. *et al.* (1985) *J. Biol. Chem.* **260**, 9976-9980Why does Coomassie® Brilliant Blue R Interact Differently with Different Proteins?