


Specification

Iodoacetamide *BioChemica*

A1666

Physical Description:	Solid
Product Code:	A1666
Product Name:	Iodoacetamide <i>BioChemica</i>
Headline Comment:	• May darken in storage!
Specifications:	Assay (HPLC): min. 99 % Water: max. 6 %
Hazard pictograms	
WGK:	1
Storage:	2-8°C
Signal Word:	Attention
GHS Symbols:	GHS07
H Phrases:	H315 H317 H319 H335
P Phrases:	P280 P302+P352 P304+P340 P305+P351+P338
Molecular Formula:	C ₂ H ₄ INO
M:	184.96 g/mol

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Specification

Iodoacetamide BioChemica

A1666

CAS:	144-48-9
EINECS:	205-630-1
CS:	29241900
<p>Comment</p> <p>Iodoacetamide and Iodoacetate are well suited for the investigation of the active center of enzymes with sulfhydryl groups in its active center. The reaction conditions for the carboxymethylation may vary with every protein, which shall be modified by the corresponding reagents. Especially the optimal pH value and the concentration of the reagent (iodoacetamide, iodoacetate, chloroacetamide, chloroacetic acid) may vary (see table II in ref. 1). In SDS polyacrylamide gel electrophoresis, the oxidation of sulfhydryl-containing proteins during electrophoresis may become a problem. Oxidation may non-reproducibly change the mobility of these proteins. Relief may come through reduction with DTE and subsequent alkylation with iodoacetamide of the proteins before loading the gel (2). Iodoacetamide is used as a protease inhibitor too. It is added to several protease inhibitor cocktails in a concentration of 1 - 5 mM. At a concentration of 4 mM, e. g. the activity of the serine protease Cerastocytin has been inhibited by more than 95 % (3). Please note, that this inhibitor, like iodoacetate and N-Ethylmaleimide, will be inactivated by reducing agents such as DTT. Especially for the preparation of plant cell extracts, high concentrations of reducing agents have to be added to the homogenizing buffers (4). Iodoacetamide may become yellowish upon longer storage!</p>	
<p>Bibliography</p> <p>(1) Gurd, F.R. (1972) <i>Methods Enzymol.</i> 25, 424-438 Carboxymethylation. (2) Lane, L.C. (1978) <i>Anal. Biochem.</i> 86, 655-664A simple method for stabilizing protein-sulfhydryl groups during SDS-gel electrophoresis. (3) Marrakchi, N. <i>et al.</i> (1995) <i>Biochim. Biophys. Acta</i> 1244, 147-156 Cerastocytin, a new thrombin-like platelet activator from the venom of the tunisian viper <i>Cerastes cerastes</i>. (4) Gegenheimer, P. (1990) <i>Methods Enzymol.</i> 182, 174-193 Preparation of extracts from plants.</p>	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340