



PRODUCT CODE: 253500

Fixing B-5 for clinical diagnosis
for microscopy

TARIC 3822 00 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,051 •

BIBLIOGRAPHY: Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore , 19 (1981) ; R. García del Moral , "Laboratoriode AnatomíaPatológica", 1ª Ed., McGraw-Hill , 617 (1993) ; E.B. Prosphet , B.Mills , J. B. Arrington , L. H. Sobin , "Laboratory Methods inHistotechnology", Armed Forces Institute of Pathology , 26 (1992) •

HAZARDOUS:



H: H330 • H310 • H300 • H373 • H411 •

P: P260 • P262 • P264 • P270 • P271 • P501 • P273 • P280 • P284 • P301+P310 • P302+P350 • P304+P340 • P310 • P314 • P320 • P321 • P322 • P330 • P361 • P363 • P391 • P403+P233 • P405 •

TRANSPORT REGULATIONS: UN: 2810 • ADR: 6.1/II • IMDG: 6.1/II • IATA: 6.1/II • PAX: 654 • CAO: 662 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~1,051 kg 1kg~0,951 l

COMPOSITION:

| | |
|--------------------------|--------|
| Mercury(II) Chloride | 6 g |
| Sodium Acetate 3-hydrate | 2,06 g |
| Water (s.q.m) | 100 ml |

SPECIFICATIONS:

Identity :
Identity

passes test

Ed.: 2 . Vig.: 30.08.2017 .

Prod.: 253500