



PRODUCT CODE: 253453

Weigert's Hematoxylin solution A for clinical diagnosis
for microscopy

TARIC 3203 00 10 00

PHYSICAL DATA: liquid, Clear, Orange, Miscible with alcohol and water • D 20/4 0,812 • Flash P.:9 °C •

INSTRUCTIONS FOR USE: Mix equal parts of solution A and solution B. Stain with this Weigert's working solution for 10 minutes.

BIBLIOGRAPHY: Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore (1981) ; H. J. Conn , "Biological Stains", 9th Ed., Williamsand Wilkins , Baltimore , 609 (1977) ; R. García del Moral,"Laboratorio de Anatomía Patológica", 1ª Ed., McGraw-Hill , 133 (1993); S. Frankel , S. Reitman , A. C. Sonnenwirth , "Gradwohl#s clinicallaboratory methods and diagnosis" , 7th Ed., The C. V. Mosby Company 1, 943 (1970) •

HAZARDOUS: LD50 oral rat 7.060 mg/kg • VLA-ED (ethanol) 1.000 ppm1.910 mg/m3



H: H225 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P280 • P303+P361+P353 • P370+P378 • P403+P235 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,812 kg 1kg~1,232 l

COMPOSITION:

Hematoxylin	1 g
Ethanol absolute	100 ml

SPECIFICATIONS:

Identity :
Identity

passes test

Ed.: 3 . Vig.: 11.12.2015 .

Prod.: 253453