



PRODUCT CODE: 253454

Weigert's Hematoxylin solution B for clinical diagnosis
for microscopy

TARIC 3203 00 10 00

PHYSICAL DATA: liquid, Clear, Greenish-yellow, Miscible with water • D 1,013 •

INSTRUCTIONS FOR USE: Mix equal parts of solution A and solution B. Stain with this Weigert's working solution for 10 minutes.

BIBLIOGRAPHY: Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore (1981) ; H. J. Conn , "Biological Stains", 9th Ed., Williamsand Wilkins , Baltimore , 609 (1977) ; R. García del Moral,"Laboratorio de Anatomía Patológica", 1ª Ed., McGraw-Hill , 133 (1993); S. Frankel , S. Reitman , A. C. Sonnenwirth , "Gradwohl#s clinicallaboratory methods and diagnosis" , 7th Ed., The C. V. Mosby Company 1, 943 (1970) •

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,013 kg 1kg~0,987 l

COMPOSITION:

Iron(III) Chloride 30% aqueous solution.	4 ml
Hydrochloric Acid 35%	1 ml
Water (s.q.m)	100 ml

SPECIFICATIONS:

Identity :	
Identity	passes test

Ed.: 3 . Vig.: 11.12.2015 .

Prod.: 253454