



**PRODUCT CODE: 256991**

**Harris Hematoxylin modified solution for clinical diagnosis**

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples., Sanitary product for in vitro diagnostic class A.

---

TARIC 3822 00 00 00

**PHYSICAL DATA:** liquid, garnet, Miscible with water

**BIBLIOGRAPHY:** Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore (1981) ; H. J. Conn , "Biological Stains", 9th Ed., Williamsand Wilkins , Baltimore , 472 (1977) ; R. García del Moral,"Laboratorio de Anatomía Patológica", 1<sup>a</sup> Ed., McGraw-Hill , 129 (1993); E.B. Prosphet , B. Mills , J. B. Arrington , L. H. Sabin , "Laboratory Methods in Histotechnology", Armed Forces Institute ofPathology , 56 (1992) ; S. Frankel , S. Reitman , A. C. Sonnenwirth , "Gradwohl's clinical I •

**HAZARDOUS:**

**WEIGHT/VOLUME INFORMATION:** 1l~1,095 kg      1kg~0,913 l

**OBSERVATIONS:** Filter any sediment present in the solution. •

**SPECIFICATIONS:**

Suitability for hematoxylin - eosin stain	passes test
---	-------------

Ed.: 3 . Vig.: 14.02.2018 .

Prod.: 256991