



PRODUCT CODE: 416116

PALCAM Listeria Selective Enrichment (Supplement) for microbiology

Physical Data: Lyophilized of orange colour · Preparation: Dissolve the contents of one vial of PALCAM Listeria Selective Supplement (Code 416116) in 2 ml sterile distilled water. Mix gently till dissolved and add, aseptically, to 500 ml of medium PALCAM Listeria Agar Base (Code 415380) autoclaved and cooled to 45-50°C. Mix well and distribute into sterile Petri dishes ·

TARIC 3821 00 00 00

PHYSICAL DATA:

BIBLIOGRAPHY: Ottaviani, F., Ottaviani, M. and Agosti, M (1987) Quimper Froid Symposium Proceedings, P6
A.D.R.I.A Quimper (F) 16-18 June ISO 11290-1:2017 Horizontal method for the detection and enumeration of Listeria monocytogenes Part 1: Detection Method. •

HAZARDOUS:



H: H318 • H412 •

P: P273 • P280 • P305+P351+P338 • P310 • P501 •

OBSERVATIONS: Storage between +2 and +8°C, in a dry place, away from direct light. •

Composition (mg/vial):

Polymixin B Sulphate	5,0
Ceftazidim	10,0
Acryflavine	2,5

Ed.: 4 . Vig.: 07.05.2018 .

Prod.: 416116