



PRODUCT CODE: 433744

Baird-Parker Agar (ISO 6888) (Contact Plate) for microbiology

Physical Data: Prepared media of opalescent yellow colour · Preparation: Get the culture media in contact with the test surface. Incubate at 37°C for 48 hours. Count the number of black colonies surrounded by a transparent halo and express the result as cfu/plate ·

TARIC 3821 00 00 00

SYNONYMS: Baird-Parker, Selective Agar, Medium O, Selective Agar for Staphylococci according to Baird-Parker •

PHYSICAL DATA:

BIBLIOGRAPHY: ISO 6888-1:1999 Microbiology of food and animal feeding stuffs: Horizontal method for enumeration of coagulase-positive staphylococci. Part1: Technique using Baird-Parker agar medium • USP 38 • Atlas R.M. & Parks L.C., Handbook of Microbiological Media ,112 (1993) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Egg Yolk Emulsion	47,0
Meat Extract	5,0
Yeast Extract	1,0
Glycine	12,0
Lithium Chloride	5,0
Potassium Tellurite	0,1
Sodium Pyruvate	10,0
Tryptone	10,0
Agar	20,5
pH: 6.9±0.2	

Ed.: 3 . Vig.: 15.01.2011 .

Prod.: 433744