

**PRODUCT CODE: 413744****Baird-Parker Agar Base (Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of lighttoasted colour · Preparation: Prepare a suspension of 63 g in 1 l of distilled water. Heat and stir to boiling point and boil for about 1 minute. Sterilize at 121°C for 15 minutes. Allow to cool to 45°C. Add 50 ml of EggYolk Tellurite Emulsion, homogenize and distribute into sterile Petri dishes ·

TARIC 3821 00 00 00

SYNONYMS: Baird-Parker, Selective Agar, Medium O, Selective Agar for Staphylococci according to Baird-Parker

PHYSICAL DATA:

BIBLIOGRAPHY: ISO 6888-1:1999 Microbiology of food and animal feeding stuffs: Horizontal method for enumeration of coagulase-positive staphylococci. Part1: Technique using Baird-Parker agar medium • USP **39** • Atlas R.M. & Parks L.C., Handbook of Microbiological Media ,112 (1993) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place. •

Composition (g/l):

Meat Extract	5,0
Yeast Extract	1,0
Glycine	12,0
Lithium Chloride	5,0
Digest Pancreatic of Casein	10,0
Sodium Pyruvate	10,0
Agar	20,0 g
pH: 6.8±0.2	

Ed.: 5 . Vig.: 15.01.2011 .

Prod.: 413744