


Specification

Hygromycin B solution

A2175

Physical Description:	Liquid
Product Code:	A2175
Product Name:	Hygromycin B solution
Specifications:	Activity: approx. 46 KU/mL Composition: Hygromycin B: approx. 41.2 mg/mL
Hazard pictograms	
UN:	2810
Class/PG:	6.1/III
ADR:	6.1/III
IMDG:	6.1/III
IATA:	6.1/III
WGK:	1
Storage:	-20°C
Signal Word:	Danger
GHS Symbols:	GHS05 GHS06 GHS08
H Phrases:	H302+H312 H317 H318

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Specification

Hygromycin B solution

A2175

	H330
	H334
P Phrases:	P280 P284 P301+P310 P305+P351+P338
Molecular Formula:	$C_{20}H_{37}N_3O_{13}$
M:	527.53 g/mol
CAS:	31282-04-9
EINECS:	250-545-5
CS:	38220000
Comment	<p>Hygromycin B belongs to the aminoglycoside antibiotics. It inhibits the bacterial translation by interacting with the 16S-rRNA. In eukaryotes the antibiotic inhibits the splicing of the group I-introns by binding to a specific RNA motive. This leads to reading errors and an inhibition of the translocation, probably by release of the peptidyl-tRNA from the acceptor site (3). The viral replication will be inhibited by antibiotics of these group too. In molecular biology, hygromycin B is used for the selection of hygromycin-resistant cells after transformation with a plasmid, that codes for the <i>E. coli</i> hygromycin-phosphotransferase (hph). A number of vectors has been developed, to enable the application in different organism (fungi): Egelhoff, T.T. <i>et al.</i> (1989) <i>Mol. Cell. Biol.</i> 9, 1965-1968; Hamada, W. <i>et al.</i> (1994) <i>Curr. Genet.</i> 26, 251-255; mammalian cells: Belt, P.B.G.M. <i>et al.</i> (1989) <i>Gene</i> 84, 407-417; Gäken, J. <i>et al.</i> (1992) <i>BioTechniques</i> 13, 32-34; Norman, J.A. <i>et al.</i> (1992) <i>Mol. Pharmacol.</i> 41, 53-59; Hubbard, S.C. <i>et al.</i> (1994) <i>J. Biol. Chem.</i> 269, 3717-3724; plants: Dale, E.C. & Ow, D.W. (1991) <i>Proc. Natl. Acad. Sci. USA</i> 88, 10558-10562; Ma, H. <i>et al.</i> (1992) <i>Gene</i> 117, 161-167). Depending on the organism and cell type, the active concentration ranges from 20 µg/ml to 150 µg/ml. Usually, the higher concentration is applied (for HeLa-cells up to 300 µg/ml). It is recommended to test for the optimal concentration. The quantity of hygromycin applied can be reduced by increasing the pH of the medium. At higher pH values, cells are more sensitive. Lower salt concentrations will help to save hygromycin too (LB medium / LB agar with 5 instead of 10 g NaCl per litre). Stability: Hygromycin B-solutions (e. g. in HEPES buffer) are stable at +4°C for at least one year. Informations for the stability at 37°C are ranging from one month to the recommendation to change the medium every third day. In agar plates the stability is approx. 4 weeks. The shipment will be at room temperature. The substance is sensitive to high concentrations of acids. Diluted acids may be tolerated for a short time. Caution: Hygromycin B is very toxic and probably a carcinogen. Do not breathe dust and avoid contact with skin. Wear suitable protective clothing, gloves and eye protection.</p>

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Specification

Hygromycin B solution

A2175

Bibliography

(1)Gritz, L. & Davies, J. (1983) *Gene* **25**, 179-188Plasmid-coded Hygromycin B-resistance\ The sequence of the Hygromycin B-Phosphotransferase gene. (2)Carrasco, L. & Vázquez, D. (1993) *Viral Translation Inhibitors in Antibiotics VI* Pages 279 (ed. F.E. Hahn) Springer-Verlag Berlin Heidelberg New York Tokyo. (3)Moazed, D. & Noller, H.F. (1987) *Nature* **327**, 389-394Interaction of antibiotics with functional sites in the 16S rRNA.

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340