

Specification

Loading buffer DNA II

A2571

Product Code:	A2571
Product Name:	Loading buffer DNA II
Short Description:	usage: 6X Gel loading buffer for DNA samples 6fold concentrated solution for agarose and acrylamide gel electrophoresis
Specifications:	Composition: Bromophenol blue sodium salt: 0.25 % Ficoll® 400: 15 % Xylene cyanol FF: 0.25 %
WGK:	1
Storage:	RT
CS:	38220000
Bibliography (1) Sambrook, J. & Russell, D.W. (2001) <i>Molecular Cloning</i> : A Laboratory Manual, 3rd Edition. Pages 1.73, 12.71, A8.18. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340