



**PRODUCT CODE: 413748**

**Brilliant Green Bile 2% Broth (ISO 4831, ISO 4832)(Dehydrated Culture Media) for microbiology**

Physical Data: Dehydrated culture media of greenish beige colour · Preparation: Dissolve 40 g in 1 l of distilled water. Distribute into tubes containing Durham bell, in 10 ml amounts for testing 1 ml or less of sample. If the sample to be analysed is bigger, prepare the medium double or triple strength. Autoclave for 15 minutes at 121°C ·

---

TARIC 3821 00 00 00

**SYNONYMS:** BGBL Broth, Brilliant Green Bile Broth, Brilliant Green Lactose Bile Broth •

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** Standard Methods for the examination of Water and Wastewater. 15th Ed. APHA. (1981). (ISO 4831, ISO 4832). •

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Ox Bile, Dehydrated	20,0 g
Brilliant Green	0,0133
Lactose	10,0
Gelatin Peptone	10,0
pH: 7.2±0.2	

