



**PRODUCT CODE: 416276**

**Orange Serum Agar (Dehydrated Culture Media) for microbiology**

Physical Data : Dehydrated culture media of beige colour · Preparation: Make a suspension of 40 g in 1 l of distilled water. Heat and stir constantly to boiling point and boil for 1 minute. Sterilize at 118°C for 15 minutes and distribute into sterile Petri dishes. Do not overheat ·

---

TARIC 3821 00 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** American Public Health Association (1976), Compendium of Methods for the Microbiological examination of Food, APHA Inc. Washington DC, Hay G. L. (1951), Proc. Florida State Hort. Soc., 94th Ann. Murdock D. I. and Brokaw C. H. (1958), Food Tech., 12, 573-576 •

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Casein Peptone	10,0
Yeast Extract	3,0
Orange Extract	5,0
Dextrose	4,0
Potassium di-Hydrogen Phosphate	3,0
Agar	15,0
pH: 5.5±0.2	

