


Specification

AppliClear (Xylene substitute)

A4632

Physical Description:	Liquid
Product Code:	A4632
Product Name:	AppliClear (Xylene substitute)
Specifications:	Density (15°C): 0.746 - 0.765 g/cm ³ n _{20°C/D} : 1.413 - 1.471 Boiling range: 140 - 165°C Aromates: max. 0.01 % Benzene: max. 0.0002 %
Hazard pictograms	
UN:	3295
Class/PG:	3/III
ADR:	3/III
IMDG:	3/III
IATA:	3/III
WGK:	1
Storage:	RT
Signal Word:	Danger

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340

Specification

AppliClear (Xylene substitute)

A4632

GHS Symbols:	GHS02 GHS07 GHS08 GHS09
H Phrases:	H225 H304 H336 H411
P Phrases:	P210 P240 P273 P301+P310 P331 P403+P235
CAS:	64742-49-0
EINECS:	265-151-9
CS:	27101221
Index Nr.:	649-328-00-1
Comment	AppliClear (xylene substitute) is used instead of xylene in histological work. It may be used as an intermediate for paraffin embedding, in the preparation for mounting, or in dewaxing of paraffin sections. AppliClear evaporates considerably faster compared to xylene. Thus, the solution is particularly suitable for automated use. When slides are operated by hand a rapid mode of operation is necessary; jars or glass chambers should be kept covered whenever possible.

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340