



PRODUCT CODE: 395465

Solution B for turbidity primary standard
for preparing formazin primary standard

TARIC 3822 00 00 00

PHYSICAL DATA: Miscible with water

INSTRUCTIONS FOR USE: Preparation of 4000 NTU Formazine standard solution: mix 5 ml of Turbidity Standard solution A with 5 ml of Turbidity Standard solution B. Allow to stand for 24 hours. Shake well before use.

BIBLIOGRAPHY: BP.2018 • Ph. Eur. 9.0 13 • Royal Decree I15 • Recueil des méthodes d'analyse et d'appréciation des vins: Turbidité des vins, Standard Methods for the Examination of water and wastewater 18th Ed., 1992, 2-9 •

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,021 kg 1kg~0,979 l

COMPOSITION:

Hexamethylenetetramine	10,0 g
Water (s.q.m)	100,0 ml

SPECIFICATIONS:

Assay (as (CH ₂) ₆ N ₄) (w/v)	9,9 - 10,1 %
--	--------------

